

ACTION[®] SUPINE POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions such as supine.

Pressure injury protection, from head to heels.

Head Positioners

► Please see our **Pediatric Product Sheet** for our extensive range of Pediatric Positioning Products.



40201

40201D

40218

Donut Head Pads

Designed to protect and cradle the patient's head, these pads afford easy access for surgical area of the head and neck. It can be used for a variety of procedures including ENT, plastic, general and ophthalmic.

40201	Adult	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40201D	Adult w/ Center Dish	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40218	Adult Wide	8 OD x 4 1/2 ID x 1 3/4"	20.5 OD x 11.4 ID x 4.5 cm



Contoured Head Pad

The pads contoured design provides optimum stability of the patient's head while in the supine position.

40213	7 x 7 x 3 1/4" (1 1/2" contoured)	18 x 18 x 8 cm
--------------	-----------------------------------	----------------

Ophthalmic Headrests

The Action[®] Ophthalmic Headrest Pad includes a contoured cutout for stable head placement. The robust edges are designed to support the surgeon's hands if needed during surgery.



40209	Large	12 x 10 x 3 3/8"	31 x 25.4 x 10 cm
40210	Medium	12 x 10 x 2 3/8"	31 x 25.4 x 7.3 cm

Supine Head Pad

The Action[®] Supine Head Pad may be used for most supine positions. After placing patient's head in the pad, the centering dish conforms to the patient's head for support and pressure prevention.



40200	9 x 10 x 2 1/4"	23 x 25 x 5.7 cm
--------------	-----------------	------------------

Head Pad with Center Dish

This moderately dished, square head pad may be used to minimize pressure on the occiput, neck or head of the patient as needed.



40207	9 x 10 x 1"	23 x 25 x 2.5 cm
--------------	-------------	------------------

Ophthalmic Cradle Headrest

Designed originally for cataract surgery, this versatile pad helps prevent risk of alopecia, especially during long surgeries.



40217	Adult	9 x 10 1/2 x 2 5/8"	23 x 26.7 x 6.7 cm
--------------	-------	---------------------	--------------------

Forehead Protector

The Forehead Protector is semi-circular with a flat bottom. It is used to protect the infant, child or adult patient's forehead from pressure exerted by oxygen tubes.



40208	2 1/2 x 1 5/8 x 5/8"	6.3 x 4 x 1.6 cm
--------------	----------------------	------------------

Arm/Elbow/Hand Positioners

Ulnar Forearm Protector

Designed to be wrapped around the arm and secured with the Hook & Loop attachment as the patient's arm is tucked to their side.



40303	18 x 13 x ¼"	46 x 33 x .6 cm
--------------	--------------	-----------------

Ulnar Nerve Protector

Designed to allow a greater view of the arm while permitting easy access by the anesthesiologist. Pad is placed over elbow and securely attached with hook & loop attachment.



40304	Small	15 ½ x 6 x ½"	40.3 x 15 x 1.3 cm
40304L	Large	19 x 8 ½ x ½"	48 x 21.6 x 1.3 cm

Hand/Wrist Support

Superior hand and wrist stabilization and protection. Perfect for Carpal Tunnel, Arterial Cannulization and Dialysis and other procedures.



40321	Hand/Wrist Support	11 x 3 ½ x 2 ¼"	28 x 9 x 5.7 cm
40315	Strap	25 x 2 x ¼"	63.5 x 5 x .635 cm
40317	Strap	36 x 2 x ½"	91.4 x 5 x 1 cm

Toboggan Pad

This pad is intended for use in conjunction with arm guard shield (toboggan). Designed to prevent patient's arm or leg from extending over the side of the table and to prevent leaning on or against the patient's extremity. *The shield is only sold with the set.*



40130	Toboggan Pad	16 ¼ x 8 ¼ x ¾"	41.3 x 21 x 1 cm
S40130*	Toboggan Pad Set	16 ¼ x 8 ¼ x ¾"	41.3 x 21 x 1 cm
	Shield size =	15 ¼ x 18 ¼"	38.7 x 46.4 cm


*Warranty covers pad only.



Contoured Armboard Pads

These pads offer excellent full-arm coverage. Contoured design cradles patient's arm, protects bicep, elbows and forearms.

40307	Small	16 x 5 ¾ x 2 ¼"	41 x 14.6 x 5.7 cm
40310	Medium	20 x 5 ¾ x 2 ¼"	51 x 14.6 x 5.7 cm
40306	Large	24 x 5 ¾ x 2 ¼"	61 x 14.6 x 5.7 cm



Armboard Pads

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas.

40300	Short	15 x 4 x ¼"	38 x 10.2 x .6 cm
40301	Long	24 x 4 x ¾"	61 x 10.2 x 1 cm
40309	Wide Short	15 x 6 ½ x ¼"	38 x 16.5 x .6 cm
40308	Wide Long	24 x 6 ½ x ¾"	61 x 16.5 x 1 cm
40315	Strap	25 x 2 x ¼"	63.5 x 5 x .635 cm
40317	Strap	36 x 2 x ½"	91.4 x 5 x 1 cm



Mayo Stand Cover

Use to prevent pressure and nerve damage when limb is placed on the Mayo Stand for positioning during orthopedic procedures.



40717	21 x 18 ½ x ¾"	53.3 x 47 x 1.6 cm
--------------	----------------	--------------------

Leg/Ankle/Heel Positioners

Dome-Shaped Positioners

40603	14 x 6 x 3"	36 x 15 x 7.6 cm
40603L	20 x 6 x 3"	50.8 x 15.2 x 7.6 cm



Semi-circular design with flat bottom configuration provides pressure protection under the knees while the patient is in a supine position. Helps to protect patient against foot drop for patients in prone position. Can be used under the metatarsal or under the foot, against the board, or in the semi-fowlers position.

Trapezoid Positioner

Designed for surgical procedures in which the patient must have clear leg separation.

40630	10 x 7 x 3½"	25.4 x 17.8 x 8.9 cm
-------	--------------	----------------------



Positioning Strap

Infection Control Design Positioning Strap. Welded seams and 4-way stretch microbe-resistant fabric make this a must have for any operating room. Radiolucent and MRI safe. Hook and loop attachments secure the strap to standard tables.

This strap is not meant to restrain patient.

40713	57 x 4 x ¼"	144.8 x 10.2 x .6 cm	Adjustable range: 28 x 43"	71.12 x 109.22 cm
40713L	90 x 4 x ¼"	228.6 x 10.2 x .6 cm	Adjustable range: 59½ x 79½"	151.13 x 201.93 cm



Heel Support

40502	6 x 3 ¼" x 1 ¼"	15 x 8.3 x 4.5 cm
-------	-----------------	-------------------



Place patient's heel in cup cut-out. Slit in back of cup secures the patient's Achilles tendon area.



Heel Protector Pad

40401	12 ½ x 7 ½ x ½"	32.7 x 20 x 1.3 cm
-------	-----------------	--------------------



Pads should be used for any procedure that is lengthy and involves pressure on the patient's heel area.



Heel Comfort Cup Set

40520 (Set of 2)	4 ¼ x 3 5/8 x ½"	10.8 x 9.2 x 1.3 cm
------------------	------------------	---------------------



Provides comfort and support for heels. May be used to complement anti-embolism stockings. Appropriate for recovery and radiology. Comes complete with polymer fastening strap. Universal sizing.

Advantages of AKTON® Polymer

- Clinically proven
- Distributes weight evenly
- Distributes heat evenly
- Will not leak or bottom out
- Will not support bacterial growth
- Imaging compatible
- Fast and easy cleaning
- Reusable--2 year warranty

- Cost Effective
- Tissue Equivalent
- Hypo-Allergenic
- Free of latex, plasticizers, silicone, DEHP



Support Surfaces & Overlays

- ▶ See our **Support Surface Brochure** for more detailed information on our full range of table pads.

Support Surfaces

Benefits you can expect from every Action® support surface:

- Clinically Engineered Infection Control Design
- Imaging Compatible
- Meets ergonomic weight limits for hospital staff
- Latex, Silicone and Plasticizer Free
- Static resistive
- Easy to clean
- CAL133 compliant
- Available in multiple thicknesses
- Coordinating armboards available
- Table Pads for all major Brands



We can fabricate replacement pads to fit your specifications.

Call customer service
800.228.7763 or 301.797.1414,
or email service@actionproducts.com

Performance Series

Our most popular table pad for high turnover Operating Rooms

- 4-Way Stretch Fabric
- Use for heavy frames
- Continuous layer of Akton® polymer maximizes pressure redistribution
- Bacteria / Microbe resistant solution in fibers of fabric

Performance Plus Series

- Premium Akton® polymer top for head to toe protection
- Best shearing qualities
- 360° multi-directional stretch
- Enhances patient stability
- Evenly distributes heat
- Evenly distributes pressure
- Repairable
- Bacteria / Microbe resistant

Angel Series

- Our most budget friendly alternative to the "basic" pad
- Upgraded viscofoam/ high density foam construction
- Use with Action® Overlays
- Comfort for non-anesthetized patient procedures



Overlays

The Action® O.R. Overlay is a versatile pad used for a wide range of procedures.

40105	Small	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm
40100	Medium	46 x 20 x ½"	117 x 50 x 1.3 cm
40101	Full Length	72 x 20 x ½"	183 x 50 x 1.3 cm
40103	Perineal Cutout	46 x 20 x ½" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
40104	3-piece Set	See A, B and C	
40104A	Head	9 x 20 x ½"	22.9 x 50.8 x 1.3 cm
40104B	Torso w/ perineal cutout	39 ½ x 20 x ½"	100.3 x 50.8 x 1.3 cm
40104C	Foot w/ perineal cutout	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm

Patient Protection for Supine Procedures including:

- | | | |
|----------------|------------------|---------------|
| Abdominoplasty | Cholecystectomy | Septoplasty |
| Blepharoplasty | Hernia Repair | Tonsillectomy |
| Breast Biopsy | Knee Arthroscopy | |
| Cardiovascular | Rhytidectomies | |

