



## Airway Management

Together we care



## Airway Management

We have been the first company in Germany specializing in the development and manufacturing of fixations and accessories for ventilated patients. Till this day it has been our intention to support the nursing staff and doctors as well as ventilated patients with a wide range of medical devices. Over the years we have created various fixations applicable for ventilated patients with tracheostoma as well as orotracheal and nasotracheal intubation.

Since every human being is an individual, so are your needs – we create your own product solution.





**BePa-Clip®**  
Tracheostomy Tube Holder

BePa-Clip® is a one-piece tracheostomy tube holder for adults and adolescents. The fixation ensures a reliable hold of the tracheal cannula and offers a high wearing comfort as well as dimension stability – even during a longer period of application.

REF	Product	Description	SU
200 444	BePa-Clip® STANDARD 1	26 cm padding   hook	60
200 443	BePa-Clip® STANDARD 2	26 cm padding   hook	60
200 447	BePa-Clip® STANDARD 3	30 cm padding   hook	60
200 301	BePa-Clip® VARIO 18/H	18 cm padding   hook	60
200 302	BePa-Clip® VARIO 18/K	18 cm padding   hook and loop	60
SA 21/H	BePa-Clip® VARIO 21/H	21 cm padding   hook	60
SA 21/K	BePa-Clip® VARIO 21/K	21 cm padding   hook and loop	60
200 501	BePa-Clip® VARIO 25/H	25 cm padding   hook	60
200 502	BePa-Clip® VARIO 25/K	25 cm padding   hook and loop	60
200 601	BePa-Clip® VARIO 30/H	30 cm padding   hook	60
200 602	BePa-Clip® VARIO 30/K	30 cm padding   hook and loop	60
200 701	BePa-Clip® VARIO 36/H	36 cm padding   hook	60
200 702	BePa-Clip® VARIO 36/K	36 cm padding   hook and loop	60



**BePa-KidsClip®**  
Tracheostomy Tube Holder

BePa-KidsClip® is a soft tracheostomy tube holder for neonates and infants. The particularly soft padding with cotton coating avoids decubitus ulcer and skin irritations.

REF	Product	Description	SU
SA 7/K	BePa-KidsClip® VARIO 7/K	7 cm padding   hook and loop	60
SA 9/K	BePa-KidsClip® VARIO 9/K	9 cm padding   hook and loop	60
200 201	BePa-KidsClip® VARIO 11/H	11 cm padding   hook	60
200 202	BePa-KidsClip® VARIO 11/K	11 cm padding   hook and loop	60
200 211	BePa-KidsClip® VARIO 13/H	13 cm padding   hook	60
200 212	BePa-KidsClip® VARIO 13/K	13 cm padding   hook and loop	60
200 221	BePa-KidsClip® VARIO 15/H	15 cm padding   hook	60
200 222	BePa-KidsClip® VARIO 15/K	15 cm padding   hook and loop	60
SA 16,5/K	BePa-KidsClip® VARIO 16,5/K	16,5 cm padding   hook and loop	60



**Duo-Clip**  
Two-Part Tracheostomy Tube Holder

Duo-Clip is a two-part tracheostomy tube holder for adults and adolescents. It is infinitely variable and hereby adaptable to any neck circumference. The fixation is available in two sizes optionally with hook or hook and loop fastener.

REF	Product	Description	SU
200 801	Duo-Clip H	32,5 cm + 10 cm padding   hook	60
200 801E	Duo-Clip H, XL	32,5 cm + 32,5 cm padding   hook	60
200 802	Duo-Clip K	32,5 cm + 10 cm padding   hook and loop	60
200 802E	Duo-Clip K, XL	32,5 cm + 32,5 cm padding   hook and loop	60



**Duo-KidsClip**  
Two-Part Tracheostomy Tube Holder

Duo-KidsClip is a two-part tracheostomy tube holder for neonates and infants. The soft padding with cotton coating optimally protects the tender skin of the small patients from decubitus ulcer. Duo-KidsClip is available in one size with hook and loop fastener.

REF	Product	Description	SU
200 805	Duo-KidsClip	19,5 cm + 9,5 cm padding   hook and loop	60



**Duo-Fix**  
Two-Part Tracheostomy Tube Holder

Duo-Fix is a two-part tracheostomy tube holder for adults and adolescents. It is infinitely variable and hereby adaptable to any neck circumference. Duo-Fix is available in one size with hook and loop fastener. Due to its shape the fixation is interesting in particular for use in homecare services.

REF	Product	Description	SU
200 803	Duo-Fix	26 cm + 16 cm padding   hook and loop	60



**Trachea-Fix**  
Tracheostomy Tube Holder

Trachea-Fix is a one-piece tracheostomy tube holder for adults and adolescents.

REF	Product	Description	SU
200 442	Trachea-Fix	Hook and loop	60



**Trachea-Fix Kids**  
Tracheostomy Tube Holder

Trachea-Fix Kids is a one-piece tracheostomy tube holder for children.

REF	Product	Description	SU
200 542	Trachea-Fix Kids	Hook and loop	60





Novo-Pad

Tracheostomy Foam Dressing

Novo-Pad is a particularly soft, skin-friendly and absorbent tracheostomy foam dressing for adults and adolescents. The dressing ensures an optimal wound milieu, prevents irritation and provides the patient with advanced tracheostomy care.

REF	Product	Description	SU
70-10	Novo-Pad	8,2 x 9,5 cm (w x l)	25
70-30	Novo-Pad 2	14 x 8 cm (w x l)	25



Trach-Basic

Tracheostomy Tube Holder, SP

Trach-Basic is a one-piece tube holder for tracheal cannulae in adults and adolescents.

REF	Product	Description	SU
100 441	Trach-Basic	Fixation length min. 17 cm   max. 46 cm	30

SP = Single Pack



Trachea-Lock Set

Four-Part Tracheostomy Tube Holder

Trachea-Lock Set is an adhesive tracheostomy tube holder for adults and adolescents. It consists of two skin-friendly plasters and two fixation bands with hook and loop fasteners. The fixation is particularly suitable for patients, where it is contraindicated to apply a tube holder around the neck (e.g. trauma, neurosurgical patients).

REF	Product	Description	SU
200 440	Trachea-Lock Set	2 plasters   2 fixation bands   hook and loop	60
200 441	Trachea-Lock fixation bands	Hook and loop	60



Sim-Plex

Fixation for Tracheal Cannulae, Bite Guards and O<sub>2</sub> Masks

Sim-Plex is a one-piece elastic band with hook fastener for fixing tracheal cannulae in adults and adolescents. Sim-Plex can also be used for the fixation of bite guards in the field of endoscopy.

REF	Produkt	Description	SU
200 446	Sim-Plex	26 – 46 cm length (max.)   hook	60



TuBo-Clip

Endotracheal Tube Holder

TuBo-Clip is a one-piece endotracheal tube holder for oral or nasal application. TuBo-Clip is available in two lengths and suitable for adults, adolescents and children. The fixation is safe and easy to handle and therefore preferred to apply for initial medical aid.

REF	Product	Description	SU
20-01	TuBo-Clip	72 cm length	100
20-05	TuBo-Clip XL	85 cm length	100



Novo-Pad V

Nonwoven Tracheostomy Dressing

Novo-Pad V is a comfortably soft, absorbent nonwoven tracheostomy dressing for adults and adolescents. The polyethylene mesh coating prevents the dressing from sticking to the skin.

REF	Product	Description	SU
70-11	Novo-Pad V, 1-ply	8,2 x 9,5 cm (w x l)	500
70-12	Novo-Pad V, 2-ply	8,2 x 9,5 cm (w x l)	250



Novo-KidsPad V

Nonwoven Tracheostomy Dressing

Novo-KidsPad V is a very soft, absorbent nonwoven tracheostomy dressing for children. The polyethylene mesh coating prevents the dressing from sticking to the skin.

REF	Product	Description	SU
70-21	Novo-KidsPad V, 1-ply	6,5 x 6,5 cm	750
70-22	Novo-KidsPad V, 2-ply	6,5 x 6,5 cm	450





### TuBo-Clip S

Endotracheal Tube Holder

TuBo-Clip S is a one-piece endotracheal tube holder for oral application in adults, adolescents and children. The adhesive hook and loop strip strengthens the safe fixation, particularly when repositioning intubated patients in the prone position and during early mobilization. TuBo-Clip S is available in two lengths.

REF	Product	Description	SU
20-02	TuBo-Clip S	72 cm length	100
20-06	TuBo-Clip S, XL	85 cm length	100
20-99	TuBo-Clip S adhesive hook and loop tape	Component for tube holder	100



### TuBo-Fix

Two-Part Endotracheal Tube Holder

TuBo-Fix is a two-part endotracheal tube holder for oral application in adults, adolescents and children. In case of dirt the middle part can be replaced separately. The neck padding might stay on the patient for a longer period.

REF	Product	Description	SU
20-10	TuBo-Fix	Middle part + neck padding, 1 pce of adhesive tape	40
20-11	TuBo-Fix middle part	26 cm length, 1 pce of adhesive tape	80
20-12	TuBo-Fix middle part XL	34 cm length, 1 pce of adhesive tape	80
20-15	TuBo-Fix neck padding	32,5 cm length	80
20-10K	TuBo-Fix 1	Neck padding, middle part, 1 pce of adhesive tape glued	40
20-10-02	TuBo-Fix 2	Neck padding, middle part, 2 pcs of adhesive tape glued	40
20-10-03	TuBo-Fix 3	Neck padding, middle part 1 of pce adhesive tape, 1 pce of adhesive tape   100 pcs of adhesive tape extra	40
20-10-TS	TuBo-Fix 4	Neck padding, middle part   40 pcs of hook and loop strip	40
20-10-05	TuBo-Fix 5	Neck padding, TuBo-Lock middle part, 1 pce of adhesive tape	40
20-11-03	TuBo-Fix middle part 3	26 cm length, 1 adhesive tape   100 pcs of adhesive tape extra	80
20-12-03	TuBo-Fix middle part 3, XL	34 cm length, 1 adhesive tape   100 pcs of adhesive tape extra	80



### TuBo-Flex

Endotracheal Tube Holder

TuBo-Flex is a latex-free and non-adhesive one-piece endotracheal tube holder for oral application in adults, adolescents and children.

REF	Product	Description	SU
20-30	TuBo-Flex	80 cm length	100




### TuBo-Grip

Two-Part Endotracheal Tube Holder

TuBo-Grip is a two-part endotracheal tube holder for adults, adolescents and children. TuBo-Grip consists of a neck padding and a middle part. The neck padding is available in two different lengths. The silicone layer of the elastic strap adheres to the tube and provides good fixing abilities even with intensely salivating patients.

REF	Product	Description	SU
20-50	TuBo-Grip	Middle part + neck padding	60
20-51	TuBo-Grip XL	Middle part + neck padding XL	60
20-52	TuBo-Grip neck padding	32 cm length	60
20-53	TuBo-Grip neck padding XL	47 cm length	60
20-55	TuBo-Grip middle part	33 cm length	60




### TuBo-Guard

Bite Guard

TuBo-Guard is a bite guard for tubes in adults, adolescents and children. This bite guard impacts no pressure upon the tongue and protects the teeth as well as the complete oropharyngeal cavity of the patient. The bite guard is directly connected to the tube and does not need much space in the patient's mouth.

REF	Product	Description	SU
20-59-7	TuBo-Guard white	Bite guard, approx. ID 7.0	20
20-59-75	TuBo-Guard blue	Bite guard, approx. ID 7.5	20
20-59-8	TuBo-Guard green	Bite guard, approx. ID 8.0 + 8.5	20
20-59-9	TuBo-Guard grey	Bite guard, approx. ID 9.0 + 9.5	20




### TuBo-Guard Set

Bite Guard with Two-Part Endotracheal Tube Holder

TuBo-Guard Set is a combination of the TuBo-Grip endotracheal tube holder and the TuBo-Guard bite guard. In case of dirt the middle part can be replaced separately, the neck padding might stay on the patient for a longer period.


REF	Product	Description	SU
20-50-7	TuBo-Guard Set white	Bite guard, approx. ID 7.0   Endotracheal Tube Holder TuBo-Grip	16
20-50-75	TuBo-Guard Set blue	Bite guard, approx. ID 7.5   Endotracheal Tube Holder TuBo-Grip	16
20-50-8	TuBo-Guard Set green	Bite guard, approx. ID 8.0 + 8.5   Endotracheal Tube Holder TuBo-Grip	16
20-50-9	TuBo-Guard Set grey	Bite guard, approx. ID 9.0 + 9.5   Endotracheal Tube Holder TuBo-Grip	16
20-51-7	TuBo-Guard Set white XL	Bite guard, approx. ID 7.0   Endotracheal Tube Holder TuBo-Grip XL	16
20-51-75	TuBo-Guard Set blue XL	Bite guard, approx. ID 7.5   Endotracheal Tube Holder TuBo-Grip XL	16
20-51-8	TuBo-Guard Set green XL	Bite guard, approx. ID 8.0 + 8.5   Endotracheal Tube Holder TuBo-Grip XL	16
20-51-9	TuBo-Guard Set grey XL	Bite guard, approx. ID 9.0 + 9.5   Endotracheal Tube Holder TuBo-Grip XL	16
20-51-X	TuBo-Guard Set Mix XL	Bite guard, 4 sets each size   TuBo-Grip XL	16




REF	Product	Description	SU
25-10	TuBo-Lock Set	2 pcs of plaster fixation   middle part   fingerlift	80
25-11	TuBo-Lock Set XL	2 pcs of plaster fixation   Middle part XL   fingerlift	80
25-11-10	TuBo-Lock Set 1, XL	2 pcs of plaster fixation   Middle part with adhesive hook and loop tape	80
25-15	TuBo-Lock middle part	28 cm length, fingerlift	80
25-16	TuBo-Lock middle part XL	36,5 cm length, fingerlift	80
25-16-10	TuBo-Lock middle part 1, XL	36,5 cm length, hook and loop tape	80
25-12	TuBo-Lock Set 2	2 pcs of plaster fixation   Middle part with silicone layer	80
25-14	TuBo-Lock middle part 2	33 cm length   with silicone layer	80
25-10-01	TuBo-Lock	Plaster fixation	80



REF	Product	Description	SU
20-60	TuBo-Cuff Monitor	Cuff Pressure Manometer	12



REF	Product	Description	SU
40020	Jack Pro® Gel plus	Oral Hygiene Gel   Dosing bottle system 100 ml	12
40030	Jack Pro® Rinsing plus	Oral Hygiene Rinsing   bottle 250 ml	6
40023	Jack Pro® Gel plus	Oral Hygiene Gel   Therapy bag system 300 ml	6
40033	Jack Pro® Rinsing plus	Oral Hygiene Rinsing   Therapy bag system 300 ml	6
41003	Set	Jack Pro® Gel plus   100 ml Jack Pro® Rinsing plus   250 ml „Oral Hygiene Knigge“	1
41004	Jack Pro® 24-36-48h SET	2x Jack Pro® Gel plus 30 ml 1x Jack Pro® Liquid plus 35 ml	20

[illegible]



With the credo **Together we care** Novo Klinik-Service GmbH develops, produces and markets innovative medical aids and medical technology solutions for the intensive care, anaesthesia and care sectors. For more than 30 years, our name has been synonymous with quality, user comfort and balanced price-performance ratio.

A decorative graphic on the left side of the page consisting of a vertical stack of overlapping hexagons in various colors: blue, red, purple, orange, and grey.

## Airway Management

## Infusion and Surgical Drain Management

## Fixation and Pressure Relief Management

## Hygiene Management

## Medical Technology



**Novo Klinik-Service GmbH**  
Technologiepark West  
Haus 15  
Zum Frenser Feld 1  
50127 Bergheim  
Deutschland | Germany

Tel. +49 2271 99404-0  
Fax +49 2271 99404-100  
[info@its-novo.com](mailto:info@its-novo.com)  
[its-novo.com](http://its-novo.com)