



Infusion and Surgical Drain Management


Together we care



Infusion and Surgical Drain Management

Our Infusion and Surgical Drain Management offers you a comprehensive spectrum of safe, comfortable, user-friendly and high quality products. Derived from practice, developed in close collaboration with intensive care specialists and nurses, they are available in different sizes and/or with variable setting options – matching the requirements of a successful therapy.

Since every human being is an individual, so are your needs – we create your own product solution.





ECMO-Fix

Fixation for Catheter Tubing Systems

ECMO-Fix serves for an easy and safe fixation of catheter tubing sets during ECMO / ECLS therapy. ECMO-Fix 1 is suitable for the fixation on the head. ECMO-Fix 2 is applied for the fixation on the leg.

REF	Product	Description	SU
100-300-01	ECMO-Fix 1	Head fixation 2,5 x 72 cm (w x l)	10
100-300-02	ECMO-Fix 2	Leg fixation 4 x 72 cm (w x l)	10



Main-Lock Set

Adhesive Tape for Catheters and Drainages

Main-Lock Set is a plaster fixation for tubing systems, catheters and wound drainages. The set consists of two components: the breathable adhesive tape with fleece surface and a fixation for tubes.

REF	Product	Description	SU
90-30-15	Main-Lock Set 10	4 x 10 cm (w x l)	40
90-30-20	Main-Lock Set 14	4 x 14 cm (w x l)	40



Main-Lock 1-3

Plaster Fixation

Main-Lock is a skin-friendly plaster with hard-wearing hook and loop straps used i.a. for the securement of infusion tubing systems, arterial or central venous catheters and feeding tubes.

REF	Product	Description	SU
90-60-81	Main-Lock 1	3,1 x 7 cm (w x l)	100
90-60-82	Main-Lock 2	4 x 9 cm (w x l)	100
90-60-83	Main-Lock 3	4 x 9 cm (w x l)	100



Main-Lock 9

Plaster Fixation

Main-Lock 9 is a skin-friendly plaster with hook and loop straps used i.a. for the securement of wound drainages, tubing systems and catheters.

REF	Product	Description	SU
90-50-09	Main-Lock 9	4 x 9 cm (w x l)	100



Nasal-Fix

Nasal Fixation for Tubes

Nasal-Fix is a skin-friendly plaster fixation for the securement of nasal tubes in adolescents and adults. Nasal-Fix can be used for both single lumen (Levin) and double lumen (Salem Sump) tubes.

REF	Product	Description	SU
100-200-01	Nasal-Fix	For all common tube sizes	80





With the credo **Together we care** Novo Klinik-Service GmbH develops, produces and markets innovative medical aids and medical technology solutions for the intensive care, anaesthesia and care sectors. For more than 30 years, our name has been synonymous with quality, user comfort and balanced price-performance ratio.

A decorative graphic on the left side of the page consisting of a vertical stack of overlapping hexagons in various colors: blue, red, purple, orange, and grey.

Airway Management

Infusion and Surgical Drain Management

Fixation and Pressure Relief Management

Hygiene Management

Medical Technology



Novo Klinik-Service GmbH
Technologiepark West
Haus 15
Zum Frenser Feld 1
50127 Bergheim
Deutschland | Germany

Tel. +49 2271 99404-0
Fax +49 2271 99404-100
info@its-novo.com
its-novo.com