

BePa-Clip® VARIO

Tracheostomy Tube Holder

BePa-Clip® VARIO is a one-piece tracheostomy tube holder for adults and adolescent patients.

It is the premium fixation on the present market.



12 |

10 |

Please control alteration of neck girth constantly to pinpoint the proper size. We respond to any kind of individual product solution, like custom-made items or special sizes. Just contact to us.

- > safe, quick and comfortable handling
- > soft – super soft padding due to cotton coating
- > latex-free – hypoallergenic, protection against decubitus ulcer
- > slightly absorptive and breathable padding
- > dimensionally stable – even during a longer period of wearing
- > economic – optimal price-performance ratio



its-novo.com

	Description	SU	REF
1	BePa-Clip® VARIO 18/H 18 cm padding hook SP	30 pcs	100 301
2	BePa-Clip® VARIO 18/H 18 cm padding hook	60 pcs	200 301
3	BePa-Clip® VARIO 18/K 18 cm padding hook and loop SP	30 pcs	100 302
4	BePa-Clip® VARIO 18/K 18 cm padding hook and loop	60 pcs	200 302
5	BePa-Clip® VARIO 21/H 21 cm padding hook	60 pcs	SA 21/H
6	BePa-Clip® VARIO 21/K 21 cm padding hook and loop	60 pcs	SA 21/K
7	BePa-Clip® VARIO 25/H 25 cm padding hook SP	30 pcs	100 501
8	BePa-Clip® VARIO 25/H 25 cm padding hook	60 pcs	200 501
9	BePa-Clip® VARIO 25/K 25 cm padding hook and loop SP	30 pcs	100 502
10	BePa-Clip® VARIO 25/K 25 cm padding hook and loop	60 pcs	200 502
11	BePa-Clip® VARIO 30/H 30 cm padding hook SP	30 pcs	100 601
12	BePa-Clip® VARIO 30/H 30 cm padding hook	60 pcs	200 601
13	BePa-Clip® VARIO 30/K 30 cm padding hook and loop SP	30 pcs	100 602
14	BePa-Clip® VARIO 30/K 30 cm padding hook and loop	60 pcs	200 602
15	BePa-Clip® VARIO 36/H 36 cm padding hook SP	30 pcs	100 701
16	BePa-Clip® VARIO 36/H 36 cm padding hook	60 pcs	200 701
17	BePa-Clip® VARIO 36/K 36 cm padding hook and loop SP	30 pcs	100 702
18	BePa-Clip® VARIO 36/K 36 cm padding hook and loop	60 pcs	200 702

SP = single packed



Fields of application:

fits on all common tracheostomy plates

INSTRUCTION FOR USE

BePa-Clip® VARIO



1 | Place the tube holder around the patient's neck.


2 | Conduct the hook fastener through the flange of the tracheostomy plate or use the hook and loop fastener and fix it.

3 | Repeat this procedure on the other side.

4 | If necessary, adjust the individual length by displacing the plastic sliders.



5 | Fixation in final position.



Hook fastener