

The new standard in lithotomy positioning

# The Purple Stirrup Catalog



Ossano Scandinavia Oy Finland  
Länsitie 349  
66400 Laihia  
[www.ossano.fi](http://www.ossano.fi)




*creators of the purple stirrup™*

All KYRA® stirrups provide unparalleled pelvic site access,  
larger boot sizing and better infection control.

### Curved rod provides more ROM and working space

- The unique patent-pending design provides 22° more lithotomy range of motion than traditional stirrups
- Great for robotic surgery
- Design protects the gas spring for longer life



### Larger boot and more padding for today's patients

- Fits patient's legs into the boots, not on top of the boots
- No exposed hard plastic edges
- *Flexible lateral fin*



### VELCRO® free straps

- Medical grade silicone straps are easy to clean and snap-in securely
- **Strong but flexible**
- Non-propagating material



Ossano Scandinavia Oy Finland  
Länsitie 349  
66400 Laihia  
www.ossano.fi

## The KYRA® Comfort™ Family of Stirrups

- Designed to improve infection control with optional Secure-Lok™ straps
- Easy to clean durable silicone straps
- **Boots sized to fit today's larger patient anatomy**
- More pressure reducing padding
- 4" more working space with wider abduction range
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°
- Lightweight

### KYRA® Comfort™ 350 Stirrups

Patient Weight Capacity 350 lbs/159 kg

Stirrup Weight 13.2 lbs/6 kg

KYRA3000 KYRA Comfort 350 Stirrups w/ traditional hook & loop straps

KYRA3100 KYRA Comfort 350 Stirrups w/ Secure-Lok™ silicone straps

### KYRA® Comfort™ 500 Stirrups

Patient Weight Capacity 500 lbs/227 kg

Stirrup Weight 14.6 lbs/6.6 kg

KYRA5000 KYRA Comfort 500 Stirrups w/ traditional hook & loop straps

KYRA5100 KYRA Comfort 500 Stirrups w/ Secure-Lok™ silicone straps

### KYRA® Comfort™ 800 Stirrups

Patient Weight Capacity 800 lbs/363 kg

Stirrup Weight 14.6 lbs/6.6 kg

KYRA8100 KYRA Comfort 800 Stirrups w/ Secure-Lok™ silicone straps



KYRA COMFORT 350



KYRA COMFORT 500



KYRA COMFORT 800



### KYRA® Single™ Stirrup Boot Covers

Infection Control for your stirrups

- Protect against cross contamination
- Fits Yellow and Purple branded stirrups
- Single use, protective barrier
- For infectious patients (known or suspected)
- For patients with open wounds, abrasion or lesions
- Decrease room turnover time

**KYRA8198** KYRA Single Boot Covers for Yellofins® Stirrups (10 pack)

**KYRA8149** KYRA Single Boot Covers for KYRA 350/500/800 Stirrups (10 pack)



creators of the purple stirrup™

# A stirrup to improve infection control.

## KYRA® Clean™ 350/500 Stirrups

350 Patient Weight Capacity 350 lbs/159 kg

500 Patient Weight Capacity 500 lbs/227 kg

- Easy to clean seamless waterproof bootpads
- No VELCRO® brand fasteners
- Easy to clean durable silicone straps
- Boots sized for today's larger patient anatomy
- More pressure reducing padding



OPTIONAL EASY TO CLEAN SEAMLESS WATERPROOF PADS



350 Stirrup Weight 13.2 lbs/5.9 kg

500 Stirrup Weight 14.4 lbs/6.5 kg

KYRA3200 KYRA Clean 350 Stirrups

KYRA5200 KYRA Clean 500 Stirrups



DURABLE SILICONE STRAPS



Ossano Scandinavia Oy Finland  
Länsitie 349  
66400 Laihia  
www.ossano.fi

## KYRA® 400 Stirrups

Patient Weight Capacity 400 lbs/181 kg

- **Non-finned boot**
- Classic PAL stirrup with updates for robotic surgery
- Easy to clean durable silicone straps
- **Boots sized to fit today's larger patient anatomy**
- More pressure reducing padding
- Increased range of motion

Stirrup Weight 12.2 lbs/5.5 kg

KYRA7100 KYRA 400 Stirrups w/ Secure-Lok™ silicone straps



KYRA 400 STIRRUP

## KYRA® Pediatrics



KYRA 250 STIRRUP

## KYRA® 250 Stirrups

Patient Weight Capacity 250 lbs/113 kg

- Easy to clean durable silicone straps
- **Boots sized to fit pediatric patients**
- More pressure reducing padding
- Lightweight

Stirrup Weight 10.2 lbs/4.6 kg

KYRA7810 KYRA 250 Stirrups w/ Secure-Lok™ silicone straps



*creators of the purple stirrup™*

# KYRA® Clamps



BEST

## Secure-Release™ Clamp

- Mounts and locks anywhere on rail
- Prevents false rail engagement
- Purple button safely releases clamp
- Autoclavable

KYRA4100	Secure-Release Clamp, US
KYRA4200	Secure-Release Clamp, EU
KYRA4300	Secure-Release Clamp, UK
KYRA4400	Secure-Release Clamp, DY
KYRA4500	Secure-Release Clamp, JP



BETTER

## Select Blade Clamp

- Mounts and locks anywhere on rail
- Prevents false rail engagement
- All metal design
- Autoclavable

KYRA4600	Select Blade Clamp, US
KYRA4610	Select Blade Clamp, EU
KYRA4620	Select Blade Clamp, UK
KYRA4630	Select Blade Clamp, DY
KYRA4640	Select Blade Clamp, JP



GOOD

## Simple Blade Clamp

- Cost effective solution
- Available for most rail sizes

KYRA4700	Simple Blade Clamp, US
KYRA4710	Simple Blade Clamp, EU/UK
KYRA4740	Simple Blade Clamp, JP

## The Purple Stirrup™ Cart

- Easy to push/pull with ergonomic handles
- Stirrups secured to prevent damage
- Clamp storage compartment
- Safely holds most stirrups
- Locking casters

KYRA1000 The Purple Stirrup Cart



Ossano Scandinavia Oy Finland  
Länsitie 349  
66400 Laihia  
www.ossano.fi



## KYRA® Armboard

- Radiolucent
- Designed to prevent fluid ingress
- Ultradurable and strong material
- Eliminates gear grinding and pinch hazards
- Patient weight capacity 500 lbs

KYRA2000	KYRA Armboard (US)
KYRA2005	KYRA Armboard (EU)
KYRA2006	KYRA Armboard (UK)
KYRA2007	KYRA Armboard (DEN)
KYRA2008	KYRA Armboard (US)



## KYRA® Comfort™ Armboard Pad

- Traditional pads with soft pressure relieving foam
- Available in various thicknesses
- Fits armboards 26" x 5.5" (66 x 14 cm)

KYRA2009	KYRA Comfort Armboard Pad 1"/2.5 cm
KYRA2010	KYRA Comfort Armboard Pad 2"/5 cm
KYRA2015	KYRA Comfort Armboard Pad 3"/7.6 cm
KYRA2016	KYRA Comfort Armboard Pad 4"/10 cm



*creators of the purple stirrup™*

# Increase Your Uptime with Kyra Medical Onsite Repair Kits

90% of stirrup repair can be done at your hospital with easy to use repair kits.



## Onsite Repair Kits Reduce Maintenance Costs

Kyra Medical provides reasonably priced repair kits that allows you to make onsite stirrup repairs 90% of the time.



## Low-Cost No-Hassle Training

Kyra Medical provides online training that will certify your biomedical/facility teams at no cost to your facility.



## No Repair Cost Surprises

Kyra Medical publishes price lists for ALL anticipated stirrup repairs<sup>1</sup>. You know **EXACTLY** what you will pay for any repair.



## Reduce Costly Equipment Downtime

Only Kyra Medical allows your biomedical team to make most repairs on site; reducing costly equipment down time.



Ossano Scandinavia Oy Finland  
Länsitie 349  
66400 Laihia  
[www.ossano.fi](http://www.ossano.fi)

## Kyra Medical Warranty

All KYRA Products are covered by a 2 (two) year warranty except for pads which carry a (1) one year warranty.

\*KYRA products may be covered by one of more issued or pending patents.

508-936-3550 | 888-611-KYRA(5972)

Contact Us Today

