



NOVO®

Tracheostomy PU foam dressing

NOVO Vinea® F/ST/Kids

Tracheostomy non-woven dressing

NOVO Vinea® N/Kids

Safe Solutions




The quick way to find the right tracheostomy dressing for you

Tracheostomy dressings Adolescents and Adults			
NOVO Vinea®			
	NOVO Vinea® F ST	NOVO Vinea® F	NOVO Vinea® N
Dimensions	90 x 90 mm	90 x 90 mm	90 x 90 mm
Material	PU foam	PU foam	1-ply non-woven
Moisture barrier	✓		
Sterile	✓		
Perforation with tear protection	✓	✓	✓
Slit	✓	✓	✓


Tracheostomy dressings Kids			
NOVO Vinea® Kids			
	NOVO Vinea® F Kids ST	NOVO Vinea® F Kids	NOVO Vinea® N Kids
Dimensions	70 x 60 mm	70 x 60 mm	70 x 60 mm
Material	PU foam	PU foam	1-ply non-woven
Moisture barrier	✓		
Sterile	✓		
Perforation with tear protection	✓	✓	✓
Slit	✓	✓	✓

The right product for your needs

NOVO tracheostomy dressings offer a safe and comfortable solution for the care of tracheostomised patients in both hospital and home settings. They are characterised by their anatomically contoured shape, light weight and highly absorbent, exceptionally soft material. At NOVO, you will find the right product for a variety of care settings.



NOVO Vinea®
dressing for adults



NOVO Vinea® Kids
dressing for children

Glossary



Non-woven



PU foam



Sterile



Paediatric product

NOVO Vinea® F ST /
NOVO Vinea® F Kids ST

Sterile Tracheostomy PU foam dressing for children and adults

Product properties

- Sterile tracheal PU foam dressing with a protective film
- Secretion-absorbent on the patient side / with a moisture barrier on the outer side
- With slit for easy change even with the cannula in place
- Optimally shaped 'Kids' dressing, adapted to the neck anatomy of infants and small children
- Highly absorbent and soft material
- Lightweight
- Suitable for standard cannulas
- Easy to use
- Latex-free, Suitable for MRI
- Sterile
- Single packed / Disposable product

Patient target group

- Infants
- (Small) Children
- Adolescents
- Adults

Material

- Polyurethane foam (PU)



REF	Article	Description	SU
N14120001	NOVO Vinea® F ST	Width 9 cm, length 9 cm, for adolescents and adults	10 pieces
N14120101	NOVO Vinea® F Kids ST	Top width 7 cm, bottom width 5 cm, length 6 cm, for infants and (small) children	10 pieces



NOVO Vinea® F / NOVO Vinea® F Kids



Tracheostomy PU foam dressing for children and adults

Product properties

- Non-sterile tracheal PU foam dressing
- With slit for easy change even with the cannula in place
- Optimally shaped 'Kids' dressing, adapted to the neck anatomy of infants and small children
- Highly absorbent and soft material
- Lightweight
- Suitable for standard cannulas
- Easy to use
- Latex-free, Suitable for MRI
- Disposable product

Patient target group

- Infants
- (Small) Children
- Adolescents
- Adults

Material

- Polyurethane foam (PU)



REF	Article	Description	SU
N14110001	NOVO Vinea® F	Width 9 cm, length 9 cm, for adolescents and adults	10 pieces
N14110101	NOVO Vinea® F Kids	Top width 7 cm, bottom width 5 cm, length 6 cm, for infants and (small) children	10 pieces

NOVO Vinea® N / NOVO Vinea® N Kids



Tracheostomy non-woven dressing for children and adults

Produktmerkmale

- 1-ply non-sterile tracheal non-woven dressing
- With slit for easy change even with the cannula in place
- Optimally shaped 'Kids' dressing, adapted to the neck anatomy of infants and small children
- Highly absorbent and soft material
- Lightweight
- Suitable for standard cannulas
- Easy to use
- Latex-free, Suitable for MRI
- Disposable product

Patientenzielgruppe

- Infants
- (Small) Children
- Adolescents
- Adults

Material

- Non-woven made of viscose, polyester and polyamide



REF	Article	Description	SU
N14010001	NOVO Vinea® N	Width 9 cm, length 9 cm, for adolescents and adults	10 pieces
N14010101	NOVO Vinea® N Kids	Top width 7 cm, bottom width 5 cm, length 6 cm, for infants and (small) children	10 pieces





Novo Klinik-Service GmbH
Technologiepark West - Haus 15
Zum Frenser Feld 1
50127 Bergheim • Germany

Tel. +49 2271 99404-0
Fax +49 2271 99404-700
info@novo-med.de
novo-med.de



scan or click