


Fixation and Pressure Relief Management



Together we care



Fixation and Pressure Relief Management

Our Fixation and Pressure Relief Management offers you a comprehensive spectrum of safe, comfortable, user-friendly and high quality products. Derived from practice, developed in close collaboration with intensive care specialists and nurses, they are available in different sizes and/or with variable setting options – matching the requirements of a successful therapy.

Since every human being is an individual, so are your needs – we create your own product solution.





Hand-Fix

Hand and Foot Fixation

Hand-Fix is a soft padded hand and foot fixation for patients of every age. It is used for self-protection of patients by preventing intended or unintended manipulation of accesses, bandages or life-sustaining technical devices.

REF	Product	Description	SU
60-01	Hand-Fix 12	3 x 12 cm (w x l) band length 65 cm	30
60-01-10	Hand-Fix 12.1	3 x 12 cm (w x l) band length 75 cm	30
60-08	Hand-Fix 17	4 x 17 cm (w x l) band length 75 cm	30
60-08-10	Hand-Fix 17.1	4 x 17 cm (w x l) band length 85 cm	30
60-50	Hand-Fix Mix	14 x Hand-Fix 12, 14 x Hand-Fix 17	1
60-61	Hand-Fix 22	7 x 22 cm (w x l) band length 150 cm	16
60-63	Hand-Fix 25	7 x 25 cm (w x l) band length 150 cm	16
60-63-01	Hand-Fix 25.1	7 x 25 cm (w x l) 100 cm sewed on	16
60-63-03	Hand-Fix 25.2	7 x 25 cm (w x l) band length 200 cm	16
60-65	Hand-Fix 30	7 x 30 cm (w x l) band length 150 cm	16
60-65-01	Hand-Fix 30.2	7 x 30 cm (w x l) band length 200 cm	16



PulsOxy-Fix



Fixation for Pulse Oximeter Sensors

PulsOxy-Fix is a softly padded fixation for pulse oximeter sensors in children. PulsOxy-Fix is offered in different widths and lengths and thus applicable regardless of sensors and age.

REF	Product	Description	SU
30-14	PulsOxy-Fix 10	1,8 x 10 cm (w x l)	120
30-10	PulsOxy-Fix 12	3 x 12 cm (w x l)	80
30-10-01	PulsOxy-Fix 12.1	3 x 12 cm (w x l) no fixation point	80
30-12	PulsOxy-Fix 17	4 x 17 cm (w x l)	60



Splint-Fix

Fixing and Positioning Splint

Splint-Fix is a softly padded fixing and positioning splint for wrist, ankle and elbow joint.

REF	Product	Description	SU
60-21	Splint-Fix 9	3-4 cm x 9 cm (w x l)	20
60-22	Splint-Fix 13	4 x 13 cm (w x l)	20
60-23	Splint-Fix 16	5 x 16 cm (w x l)	20
60-11	Splint-Fix 22	9 x 22 cm (w x l)	14
60-13	Splint-Fix 33	9 x 33 cm (w x l)	10



Oxy-Care

Cushioning for Tubing Systems

Oxy-Care is a comfortable and easy cushioning for tubes such as nasal cannulae.

REF	Product	Description	SU
40-10	Oxy-Care	Ø 1,0 cm padding length 7 cm	120




Ute-Fix

Cushioning for Tubing Systems

Ute-Fix is a soft, atraumatic cushioning to avoid skin damages caused by tubes and drainages.

REF	Product	Description	SU
40-06	Ute-Fix 1	Ø 1,5 cm padding length 2 cm	60
40-05	Ute-Fix 2	Ø 2 cm padding length 2 cm	60
40-01	Ute-Fix 3	Ø 2,5 cm padding length 2 cm	60



Eye-Care



Phototherapy Eye Protector

Eye-Care is a particularly soft eye protection for phototherapy in neonates. The comfortable materials are latex-free and non-adhesive as well as gentle to the tender skin by the cotton lining on the inside of the eye protector.

REF	Product	Description	SU
60-31	Eye-Care 32	32 cm band length	15
60-32	Eye-Care 38	38 cm band length	15



Clip-Fix

Fixation for Cables and Tubes

Clip-Fix is a hook and loop strap with a metal holder, which serves for the sorting of tubing and cable systems at the bed sheet, patient's clothing or DIN rails. Clip-Fix is available in different versions and lengths.

REF	Product	Description	SU
30-30-01	Clip-Fix 1	2,5 x 11 cm (w x l)	50
30-30-02	Clip-Fix 2	2,5 x 17 cm (w x l)	50
30-30-03	Clip-Fix 3	2,5 x 23,5 cm (w x l)	40
30-30-04	Clip-Fix 1 Duo	2,5 x 12,5 cm (w x l)	50
30-30-02D	Clip-Fix 2 DIN	2,5 x 18,5 cm (w x l)	25
30-30-03D	Clip-Fix 3 DIN	2,5 x 25 cm (w x l)	20
30-30-08	Clip-Fix 2 Trio	2,5 x 17 cm (w x l)	50



Uni-Fix

Fixation for Cables and Tubes

Uni-Fix is used in all hospital areas where cables and tubes need to be fixed and sorted remotely from the patient. The one-piece fleecy band with hook and loop fastener and rectangular ring allows a multiple, safe application.

REF	Product	Description	SU
30-50-01	Uni-Fix 1	2 x 24 cm (w x l)	20
30-50-02	Uni-Fix 2	2 x 40 cm (w x l)	20



Soft-Fix

Multi-Fixation

Soft-Fix is a non-adhesive, multifunctional fixation suitable for different fields of application. For instance, it serves for the securement of infusion cannulae, probes, drainages, catheters or cables of medical devices.

REF	Product	Description	SU
30-01	Soft-Fix S	2 cm x 5 m (w x l)	12
30-02	Soft-Fix M	4 cm x 5 m (w x l)	6
30-03	Soft-Fix L	6 cm x 5 m (w x l)	4
30-04	Soft-Fix hook and loop roll	2 cm x 15 m (w x l)	1
30-05	Soft-Fix Mix	5 x S, 2 x M, 1 x L, 1 x hook and loop roll	1




Transducer-Fix

Fixation for invasive Blood Pressure Monitoring

Transducer-Fix is a soft, extra wide cuff serving for the safe fixation of transducers in adults and adolescents. The continuous and faultless invasive measurement of blood pressure is ensured with Transducer-Fix, especially when it comes to repositioning or patient transportation.

REF	Product	Description	SU
30-45	Transducer-Fix H	Horizontal hook and loop strap	36
30-45-03	Transducer-Fix V	Vertical hook and loop strap	36





With the credo **Together we care** Novo Klinik-Service GmbH develops, produces and markets innovative medical aids and medical technology solutions for the intensive care, anaesthesia and care sectors. For more than 30 years, our name has been synonymous with quality, user comfort and balanced price-performance ratio.

A vertical column of colorful geometric shapes, including hexagons and triangles, in shades of blue, red, purple, orange, and grey, positioned on the left side of the page.

Airway Management

Infusion and Surgical Drain Management

Fixation and Pressure Relief Management

Hygiene Management

Medical Technology



Novo Klinik-Service GmbH
Technologiepark West
Haus 15
Zum Frenser Feld 1
50127 Bergheim
Deutschland | Germany

Tel. +49 2271 99404-0
Fax +49 2271 99404-100
info@its-novo.com
its-novo.com