NOVO Update

NOVO Adoxa

Extremity fixation on the operating table

NEW

Product properties

- Arm fixation on the operating table
- Quick and easy to use
- Fine foam structure on the skin side provides a good hold
- Skin-friendly materials
- Welded hook and loop fastener
- Suitable for MRI
- Disposable product




Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam
- Latex-free



REF	Article	Description	SU
N60501060	NOVO Adoxa 60	Arm fixation, 60 cm long, 4 cm wide, hook and loop welded on, hook and loop length 8.5 cm	60 pieces

Novo Klinik-Service GmbH
Technologiepark West – Haus 15
Zum Frenser Feld 1
50127 Bergheim • Germany

Tel. +49 2271 99404-0
Fax +49 2271 99404-700
info@novo-med.de
novo-med.de

Rev.01 Novo-update-Adoxa_EN_202404