NOVO Update

NOVO Malva® / NOVO Malva® Kids

Underpadding and fixation

NEW

Product properties

- Underpadding and fixation of arterial lines and cross-slotting for three-way stopcock (NOVO Malva® H)
- For padding and fixation of infusions, with thumb cut-out (NOVO Malva® H Kids)
- Quick and easy to use
- Fine foam structure on the skin side provides a good grip
- Skin-friendly materials
- Welded hook and loop fastener
- Suitable for MRI
- Disposable product

Patient target group

- Children, Adolescents, Adults

Material


- Soft, laminated foam
- Latex-free




NOVO Malva® H

NOVO Malva® H Kids

Application examples



The cross slot punching of the **NOVO Malva® H** allows padded fixation of the three-way tap with very good accessibility at the same time



The slot punch for the thumb to pass through **NOVO Malva® H Kids** allows the padding to be positioned securely

REF	Article	Description		VE
N42001120	NOVO Malva® H Kids	Hand padding with slot punching Length 20 cm, width 4 cm, for children	–	60 pieces
N42001040	NOVO Malva® H	Underpadding with cross-slot punching, length 40 cm, width 6 cm for adolescents and adults	–	60 pieces

Novo Klinik-Service GmbH
Technologiepark West – Haus 15
Zum Frenser Feld 1
50127 Bergheim • Germany

Tel. +49 2271 99404-0
Fax +49 2271 99404-700
info@novo-med.de
novo-med.de

Rev.01 Novo-Malva_EN_202405